

## Appendix 2

### Molecular Biology Masterclass: One day practical

#### Materials required (per group)

3 $\mu$ g Lambda DNA (e.g. 12  $\mu$ L of 0.3 $\mu$ g/ $\mu$ L)

4 Units *Eco*R1 (1 $\mu$ L)

10X *Eco*R1 buffer (5 $\mu$ L)

4 Units *Hind*III (1 $\mu$ L)

10X *Hind*III buffer (5 $\mu$ L)

20 $\mu$ L Pre-digested  $\lambda$ /*Eco*R1 (diluted from stock so 20 $\mu$ L contains 1 $\mu$ g DNA)

20 $\mu$ L Pre-digested  $\lambda$ /*Hind*III (diluted from stock so 20 $\mu$ L contains 1 $\mu$ g DNA)

100 $\mu$ L 6X gel loading buffer

Tubes labelled (but containing only SDW)

5 $\mu$ L pTZ19u

5 $\mu$ L Control

500 $\mu$ L Competent cells – HB101 in LB overnight 2 x 50mL @ 37°C

Ice bucket

Finn pipettes 0.1–10 $\mu$ L

5–40 $\mu$ L

40–200 $\mu$ L

100–1000 $\mu$ L

Box of sterile yellow and blue tips

10mL Sterile water

10mL Luria Broth

2x Agar plates

10 Sterile Eppendorf tubes

marker pen

2 foam Eppendorf racks

Bunsen burner

Pot of Ethanol

Spreader

Demonstration transformation plates – 1 transformed plate and 1 control plate (see methods)

## **Materials required (per class)**

Water bath set at 37°C  
Water bath set at 42°C  
Chloros pots  
Mini gel rig and cast gel (see methods)  
Power pack  
Transilluminator  
S M L Gloves

## ***Molecular Biology Masterclass — Recipes***

### **TE Buffer**

10mM Tris (pH 7.5) 1.21g/L  
1mM EDTA 0.373g/L

### **TBE Buffer**

#### **Working solution**

0.089M Tris base; 0.089M boric acid; 0.002M EDTA

#### **Concentrated stock solution (5X) – use at 1X**

Tris base	54g/L	<b><u>1X</u></b>	Tris base	10.8g/L
Boric acid	27.5g/L		Boric acid	5.5g/L
EDTA (pH 8.0)	3.7g		EDTA	.74g/L

### **6X gel loading buffer**

0.25% xylene-cyanol FF  
0.25% bromophenol blue  
40% (w/v) sucrose in distilled water

### **LB Medium (Luria broth)**

Tryptone 10g/l  
Yeast extract 5g/l  
NaCl 10g/l

Adjust to pH 7.5 with NaOH; for agar add 1.5% agar No. 2

### **Agar plates for blue colonies**

For agar use LB medium containing 1.5% agar No.2.

Per 300ml of agar (i.e. 1x 500mL medical flat) add:

- (1) 3mL of 10mg/mL ampicillin in distilled water and filter sterilise via addition with a syringe.
- (2) 0.5mL of 25mg/mL X-Gal in DMF. (Check compatibility of any plastic used.) This is self-sterile.
- (3) 1mL of 8mg/mL IPTG dissolved in distilled water and filter sterilise via addition with a syringe.

## **DNA Electrophoresis**

### **Casting a gel:**

- (1) Cover the open ends of the gel tray with two layers of autoclave tape.
- (2) Insert the well-forming comb.
- (3) Add 0.7% agarose to 1xTBE buffer. 300mL TBE for a large rig and 50ml for a mini rig.)
- (4) Melt the agarose in the microwave on full power.
- (5) When the agarose is hand hot add 0.5µg/mL ethidium bromide.
- (6) Pour the agarose into the gel tray to approximately halfway up the comb.
- (7) When the gel is fully set remove the comb and autoclave tape.
- (8) Place the tray in the gel rig and cover with 1xTBE buffer. Make sure the wells are at the negative end of the rig as the DNA runs to the positive side.

### **Loading the gel:**

- (1) Carefully load 20µL of sample into the wells.
- (2) Place the lid on the rig and switch power pack on.
- (3) Run the gel rig until the samples are at least  $\frac{3}{4}$  along the gel.  
[Mini rig @ 75volts or Large rig @ 200volts]